

Mölnlycke®

Mepilex® Border 

WITH SAFETAC® TECHNOLOGY

STERILE EO

UK  
CA 0086

CE 2797


Manufacturer / Fabricant

 Mölnlycke Health Care AB  
Gamlestadvägen 3C, Box 13080  
SE-402 52 Göteborg, Sweden


[www.molnlycke.com](http://www.molnlycke.com)

## Mepilex® Border Ag

## Mepilex® Border Sacrum Ag

[www.molnlycke.com/symbols](http://www.molnlycke.com/symbols)

Single use  
Usage unique

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# Mepilex® Border Ag

Self-adherent antimicrobial soft silicone foam dressing

## Product description

Mepilex Border Ag consists of a Safetac® soft silicone wound contact layer, an absorbent polyurethane foam pad containing silver sulphate and coloured by activated carbon, a layer with super absorbent fibres, a non-woven and a breathable and water proof film.

Safetac® is a proprietary adhesive technology developed by Mölnlycke that minimises pain to patients and trauma to intact skin or wounds.

Mepilex Border Ag is a soft silicone foam dressing that absorbs exudate, maintains a moist wound healing environment and has antimicrobial properties.

Mepilex Border Ag contains silver sulphate to create an effective bacterial barrier and inactivate a wide range of wound related pathogens (bacteria and fungi), as shown *in vitro*. By reducing the number of micro-organisms, Mepilex Border Ag may also reduce odour.

Mepilex Border Ag has been shown to inactivate wound related pathogens up to 7 days *in vitro*.

## Dressing material content:

Polyurethane, silicone, polyacrylate, cotton, viscose, polyester, polyolefin, silver sulphate in the wound pad which equals to 1.2 mg/cm<sup>2</sup> silver

## Indications for use

Mepilex Border Ag is designed for the management of medium to high exuding leg and foot ulcers, pressure ulcers, malignant wounds, partial thickness burns, traumatic and surgical wounds where a moist environment, exudate handling, gentle fixation and an antimicrobial action is indicated.

Mepilex Border Ag may be used on infected wounds as part of a treatment regimen under supervision of a qualified health care professional.

Mepilex Border Ag can be used under compression bandaging.

**Precaution(s)**

- Do not use Mepilex Border Ag on patient and/or user with known hypersensitivity to the ingoing materials/components of the product.
- Clinicians/healthcare professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Frequent or prolonged use of this product may result in transient or permanent discoloration of skin.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.

**Instructions for use**

Mepilex Border Ag should be used by or under the supervision of a qualified health care professional.

Note that local hygiene procedures should be followed prior to and following the dressing change.

1. Cleanse the wound with saline solution or water according to standard clinical practice.
2. Dry the surrounding skin thoroughly.
3. Select appropriate dressing size. For best clinical outcomes, ensure that the absorbent pad covers the wound bed and overlaps the wound edges.
4. Remove the release films and apply the adherent side to the wound. Do not stretch.

Each Mepilex Border Ag dressing may be left in place for up to 7 days depending on the patient, condition of the wound and surrounding skin, or as indicated by accepted clinical practice.

Initially it might be necessary to change Mepilex Border Ag more frequently. This is due to a change in treatment regime, which can often result in an initial increase in exudation.

Mepilex Border Ag is intended for short-term use up to 4 weeks. For continued use, a re-assessment by a physician is recommended. In case of clinical infection Mepilex Border Ag does not replace the need for systemic or other adequate infection treatment.

The interaction of Mepilex Border Ag together with cleansing agents, oxidising agents such as hypochlorite solutions/hydrogen peroxide and topical treatments has not been demonstrated.

Prior to commencing radiation therapy, remove Mepilex Border Ag if product is present in the treatment area. A new dressing can be applied following treatment.

Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).

### **Special storage conditions and handling conditions**

Mepilex Border Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight.

Disposal should be handled according to local environmental procedures.

### **Other information**

For EU-customers, further information on the safety and clinical performance of Mepilex Border Ag please refer to the 'Summary of Safety and Clinical Performance' for Mepilex Border Ag. The document is published in the European Database on Medical Devices, EUDAMED, <https://ec.europa.eu/tools/eudamed>.

If any serious incident has occurred in relation to the use of Mepilex Border Ag, you should report to Mölnlycke Health Care and to your local competent authority.

Mepilex and Safetac technology are registered trademarks of Mölnlycke Health Care AB.

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